

RBM17 Antibody

Catalog No: #47187

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RBM17 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total RBM17 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human RBM17
Target Name	RBM17
Other Names	SPF45
Accession No.	Swiss-Prot#:Q96I25NCBI Gene ID:84991Gene Accssion:BC007871
Uniprot	Q96I25
GeneID	84991;
Concentration	0.8mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20C


Application Details

Immunofluorescence:1: 40-200

Images



The image is immunohistochemistry of paraffin-embedded Human lung cancer tissue using 47187(RBM17 Antibody) at dilution 1/40. (Original magnification: 200)



The image is immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using 47187(RBM17 Antibody) at dilution 1/40. (Original magnification: ?00)

Background

This gene encodes an RNA binding protein. The encoded protein is part of the spliceosome complex and functions in the second catalytic step of mRNA splicing. Alternatively spliced transcript variants have been described. Related pseudogenes exist on chromosomes 9 and 15.

Note: This product is for in vitro research use only