

SLC22A24 Antibody

Catalog No: #47205



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	SLC22A24 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total SLC22A24 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human SLC22A24
Target Name	SLC22A24
Other Names	NET46
Accession No.	Swiss-Prot#:Q8N4F4NCBI Gene ID:283238Gene Accssion:BC034394
Uniprot	Q8N4F4
GeneID	283238;
Concentration	1.1mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN3, 40% Glycerol.
Storage	Store at -20C


Application Details

Immunofluorescence:1: 25-100

Images



The image is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using 47205(SLC22A24 Antibody) at dilution 1/30. (Original magnification: ?00)



The image is immunohistochemistry of paraffin-embedded Human prostate cancer tissue using 47205(SLC22A24 Antibody) at dilution 1/30. (Original magnification: ?00)

Background

SLC22A24 belongs to a large family of transmembrane proteins that function as uniporters, symporters, and antiporters to transport organic ions across cell membranes (Jacobsson et al., 2007 [PubMed 17714910]).

Note: This product is for in vitro research use only