

TMEM248 Antibody

Catalog No: #47220

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TMEM248 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	WB, IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total TMEM248 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human TMEM248
Target Name	TMEM248
Other Names	C7orf42
Accession No.	Swiss-Prot#:Q9NWD8NCBI Gene ID:55069Gene Accssion:BC008675
Uniprot	Q9NWD8
GeneID	55069;
Calculated MW	35 kDa
Concentration	0.8mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20C


Application Details

Western blotting:1:500-2000Immunofluorescence:1: 25-100

Images



The image is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using 47220(TMEM248 Antibody) at dilution 1/25. (Original magnification: ?00)



Gel: 8%SDS-PAGE

Lysate: 40 µg, Lane 1-2: A172 and HEPG2 cell lysates

Primary antibody: TMEM248 Antibody at dilution 1/300

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 5 seconds

Background

The function of this protein remains unknown.

Note: This product is for in vitro research use only