

TMEM45A Antibody

Catalog No: #47221

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TMEM45A Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total TMEM45A protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide of human TMEM45A
Target Name	TMEM45A
Other Names	DERP7; DNAPTP4
Accession No.	Swiss-Prot#:Q9NWC5 NCBI Gene ID:55076Gene Accssion:NP_060474
Uniprot	Q9NWC5
GeneID	55076;
Concentration	0.5mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20C


Application Details

Immunofluorescence:1: 25-100

Images



The image is immunohistochemistry of paraffin-embedded Human cervical cancer tissue using 47221(TMEM45A Antibody) at dilution 1/20. (Original magnification: ?00)



The image is immunohistochemistry of paraffin-embedded Human tonsil tissue using 47221(TMEM45A Antibody) at dilution 1/20. (Original magnification: 200)

Background

High expression and prognostic value of TMEM45A has been recently reported in various types of human tumors. High levels of TMEM45A expression in tumors may be indicative of potential resistance to cancer therapy, making TMEM45A an interesting biomarker for resistance.

Note: This product is for in vitro research use only