

## TTC4 Antibody

Catalog No: #47229

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	TTC4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total TTC4 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human TTC4
Target Name	TTC4
Accession No.	Swiss-Prot#:O95801NCBI Gene ID:7268Gene Accssion:BC001276
Uniprot	O95801
GeneID	7268;
Concentration	1.6mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20C

## Application Details

Immunofluorescence:1: 40-200

## Images



The image is immunohistochemistry of paraffin-embedded Human liver cancer tissue using 47229(TTC4 Antibody) at dilution 1/45. (Original magnification: ?00)

## Background

This gene encodes a protein that contains tetratricopeptide (TPR) repeats, which often mediate protein-protein interactions and chaperone activity. The encoded protein interacts with heat shock proteins 70 and 90. Alternative splicing results in multiple transcript variants. Naturally-occurring readthrough transcription occurs from upstream gene MROH (maestro heat-like repeat family member 7) to this gene.

---

Note: This product is for in vitro research use only