

UQCR10 Antibody

Catalog No: #47233

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	UQCR10 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total UQCR10 protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide of human UQCR10
Target Name	UQCR10
Other Names	QCR9; UCRC; HSPC051; HSPC119; HSPC151; UCCR7.2
Accession No.	Swiss-Prot#:Q9UDW1 NCBI Gene ID:29796Gene Accssion:NP_037519
Uniprot	Q9UDW1
GeneID	29796;
Concentration	0.8mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20C


Application Details

Immunofluorescence:1: 20-100

Images



The image is immunohistochemistry of paraffin-embedded Human gastric cancer tissue using 47233(UQCR10 Antibody) at dilution 1/30. (Original magnification: ?00)



The image is immunohistochemistry of paraffin-embedded Human liver cancer tissue using 47233(UQCR10 Antibody) at dilution 1/30. (Original magnification: 200)

Background

UCRC is a subunit of mitochondrial complex III (ubiquinol-cytochrome c reductase; EC 1.10.2.2), which forms the middle segment of the respiratory chain of the inner mitochondrial membrane (Schagger et al., 1995 [PubMed 8592474]).

Note: This product is for in vitro research use only