

RRAS2 Antibody

Catalog No: #40328



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	RRAS2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total RRAS2 protein.
Immunogen Type	Protein
Immunogen Description	Full length fusion protein
Target Name	RRAS2
Other Names	TC21
Accession No.	Swiss-Prot:P62070Gene Accssion:BC013106
Uniprot	P62070
GeneID	22800;
SDS-PAGE MW	23KD
Concentration	0.9mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN3, 40% Glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500-1:2000

Images



Gel: 12%SDS-PAGE
Lysates (from left to right): NIH/3T3 and A431 cell
Amount of lysate: 40ug per lane
Primary antibody: 1/200 dilution
Secondary antibody dilution: 1/8000
Exposure time: 15 seconds

Background

This gene encodes a member of the R-Ras subfamily of Ras-like small GTPases. The encoded protein associates with the plasma membrane and may function as a signal transducer. This protein may play an important role in activating signal transduction pathways that control cell proliferation.

Mutations in this gene are associated with the growth of certain tumors.

Note: This product is for in vitro research use only