

S100PBP Antibody

Catalog No: #40329

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	S100PBP Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total S100PBP protein.
Immunogen Type	Protein
Immunogen Description	Full length fusion protein
Target Name	S100PBP
Other Names	S100PBPR
Accession No.	Swiss-Prot:Q96BU1Gene Accssion:BC015175
Uniprot	Q96BU1
GeneID	64766;
SDS-PAGE MW	46KD
Concentration	1.4mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:200-1:1000

Images



Gel: 8%SDS-PAGE
Lysate: 40ug Human fetal brain tissuePrimary antibody: 1/350
dilution
Secondary antibody dilution: 1/8000
Exposure time: 1 minute

Background

This gene encodes a protein that was originally identified by its interaction with S100 calcium-binding protein P. Expression of this protein has been reported to be associated with pancreatic ductal adenocarcinoma. Alternatively spliced transcript variants have been described.?

Note: This product is for in vitro research use only