

ARL4 antibody

Catalog No: #22527

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ARL4 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human ARL4A
Target Name	ARL4
Accession No.	NCBI Gene ID: 10124NCBI mRNA#: NM_001037164NCBI Protein#: NP_001032241
Uniprot	P40617
GeneID	10124;
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 23kd

Western blotting: 1:500-1:3000

Immunohistochemistry: 1:100-1:250

Images



Western blot analysis ARL4A using HeLa whole cell lysates

Background

ADP-ribosylation factor-like 4A is a member of the ADP-ribosylation factor family of GTP-binding proteins. ARL4A is similar to ARL4C and ARL4D and each has a nuclear localization signal and an unusually high guanine nucleotide exchange rate. ARL4A is located in both the nuclear and extranuclear cell compartments. Multiple transcript variants encoding the same protein have been found for this gene. [provided by RefSeq]

Note: This product is for in vitro research use only