

KIF1B (Phospho-Ser1057) Antibody

Catalog No: #12992

Package Size: #12992-1 50ul #12992-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	KIF1B (Phospho-Ser1057) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody is from purified rabbit serum by affinity purification via sequential chromatography on phospho- and non-phospho-peptide affinity columns.
Applications	WB
Species Reactivity	Hu, Ms
Specificity	Phospho-KIF1B (Phospho-Ser1057) Antibody detects endogenous levels of KIF1B only when phosphorylated at Ser1057
Immunogen Type	peptide
Immunogen Description	A synthesized peptide derived from human KIF1B around the phosphorylation site of Ser1057
Target Name	KIF1B
Other Names	Charcot Marie Tooth neuropathy 2A; CMT 2; CMT 2A; CMT 2A1; CMT2 A; CMT2 A1; CMT2; CMT2A 1; CMT2A; CMT2A1; D4Mil1e; HMSN II; HMSNII; HMSNII, hereditary motor sensory neuropathy II; KIF 1B; KIF1 B; KIF1B p130; KIF1B p204; KIF1Bp130; KIF1Bp204; kinesin family member 1B; Kinesin like protein KIF1B; Klp;
Accession No.	Swiss-Prot#:O60333NCBI Gene ID:23095
Uniprot	O60333
GeneID	23095;
Calculated MW	204
Concentration	1mg/ml
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20 °C.Stable for 12 months from date of receipt

Application Details

WB: 1:500-1:2000

Images



Western blot analysis of extracts of Mouse brain? tissue sample,using Phospho-KIF1B (Ser1057) Antibody.(12992)

Note: This product is for in vitro research use only