

CRYAB (Phospho-Ser19) Antibody

Catalog No: #13039



Package Size: #13039-1 50ul #13039-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CRYAB (Phospho-Ser19) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody is from purified rabbit serum by affinity purification via sequential chromatography on phospho- and non-phospho-peptide affinity columns.
Applications	WB
Species Reactivity	Hu, Ms, Rt
Specificity	Phospho-CRYAB (Phospho-Ser19) Antibody detects endogenous levels of CRYAB only when phosphorylated at Ser19
Immunogen Type	peptide
Immunogen Description	A synthesized peptide derived from human CRYAB around the phosphorylation site of Ser19
Target Name	CRYAB
Other Names	AACRYA; Alpha B crystallin; Alpha crystallin B chain; Alpha(B)-crystallin; Alpha-crystallin B chain; CRYA2; Cryab; CRYAB_HUMAN; Crystallin alpha B; Crystallin alpha polypeptide 2; CTPP2; Heat shock 20 kD like protein; Heat shock protein beta 5; Heat shock protein beta-5; HspB5; Renal carcinoma antigen NY REN 27; Renal carcinoma antigen NY-REN-27; Rosenthal fiber component;
Accession No.	Swiss-Prot#:P02511NCBI Gene ID:1410
Uniprot	P02511
GeneID	1410;
Calculated MW	20
Concentration	1mg/ml
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20 °C.Stable for 12 months from date of receipt

Application Details

WB: 1:500-1:2000

Images



Western blot analysis of Phospho-CRYAB (Ser19) in lysates of HepG2 cells treated with nocodazole 1 µg/ml 16h, using Phospho-CRYAB (Ser19) Antibody(13039).

Note: This product is for in vitro research use only