

ATP1 alpha1/Na+K+ ATPase1 (Phospho-Tyr542) Antibody

Catalog No: #13060



Package Size: #13060-1 50ul #13060-2 100ul

Orders: order@signalwayantibody.com


Support: tech@signalwayantibody.com

Description

Product Name	ATP1 alpha1/Na+K+ ATPase1 (Phospho-Tyr542) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody is from purified rabbit serum by affinity purification via sequential chromatography on phospho- and non-phospho-peptide affinity columns.
Applications	WB
Species Reactivity	Hu, Ms, Rt
Specificity	Phospho-ATP1 alpha1/Na+K+ ATPase1 (Phospho-Tyr542) Antibody detects endogenous levels of ATP1 alpha1/Na+K+ ATPase1 only when phosphorylated at Tyr542
Immunogen Type	peptide
Immunogen Description	A synthesized peptide derived from human ATP1 alpha1/Na+K+ ATPase1 around the phosphorylation site of Tyr542
Target Name	ATP1 alpha1/Na+K+ ATPase1
Other Names	ATPase Na+/K+ transporting alpha; adenosinetriphosphatase; AT1A1_HUMAN; ATP1A1; ATP1A4; ATP1AL2; ATP1B; ATP1B1; ATPase Na+/K+ transporting alpha 1 polypeptide; ATPase Na+/K+ transporting alpha 4 polypeptide; ATPase Na+/K+ transporting beta 1 polypeptide; ATPase, Na+/K+ transporting, alpha polypeptide-like 2; ATPase, Na+/K+ transporting, beta 1 polypeptide; Beta 1-subunit of Na(+),K(+)-ATPase; Na(+)/K(+) ATPase alpha-1 subunit; Na(+)/K(+) ATPase alpha-4 subunit; Na+, K+ ATPase alpha subunit; Na+/K+ ATPase 1; Na+/K+ ATPase 4; Na+/K+ ATPase, alpha-D polypeptide; Na, K-ATPase beta-1 polypeptide; Na, K-ATPase, alpha-A catalytic polypeptide; Na,K-ATPase catalytic subunit alpha-A protein; Na,K-ATPase subunit alpha-C; polypeptide-like 2; Sodium pump 1; sodium pump 4; Sodium pump subunit alpha-1; sodium pump subunit alpha-4; sodium pump subunit beta-1; sodium-potassium ATPase catalytic subunit alpha-1; sodium-potassium ATPase catalytic subunit alpha-4; sodium-potassium ATPase subunit beta 1 (non-catalytic); sodium-potassium ATPase, alpha 4 polypeptide; sodium-potassium-ATPase, alpha 1 polypeptide; Sodium/potassium transporting ATPase alpha 1 chain; Sodium/potassium transporting ATPase subunit beta 1; sodium/potassium-dependent ATPase beta-1 subunit; sodium/potassium-transporting ATPase alpha-4 chain; sodium/potassium-transporting ATPase beta-1 chain; Sodium/potassium-transporting ATPase subunit alpha-1; sodium/potassium-transporting ATPase subunit alpha-4;
Accession No.	Swiss-Prot#:P06685NCBI Gene ID:24211
Uniprot	P06685
GeneID	24211;
Calculated MW	113
Concentration	1mg/ml
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20 °C.Stable for 12 months from date of receipt

Application Details

Images



Note: This product is for in vitro research use only