

EPRS (Phospho-Tyr377) Antibody

Catalog No: #13065



Package Size: #13065-1 50ul #13065-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	EPRS (Phospho-Tyr377) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody is from purified rabbit serum by affinity purification via sequential chromatography on phospho- and non-phospho-peptide affinity columns.
Applications	WB
Species Reactivity	Hu, Ms
Specificity	Phospho-EPRS (Phospho-Tyr377) Antibody detects endogenous levels of EPRS only when phosphorylated at Tyr377
Immunogen Type	peptide
Immunogen Description	A synthesized peptide derived from human EPRS around the phosphorylation site of Tyr377
Target Name	EPRS
Other Names	Bifunctional aminoacyl tRNA synthetase; Bifunctional aminoacyl-tRNA synthetase; Bifunctional glutamate/proline tRNA ligase; Cell proliferation-inducing gene 32 protein; DKFZp313B047; EARS; Eprs; GLNS; Glu pro tRNA synthetase; GLUPRORS; GluRS; Glutamate tRNA ligase; Glutametyl prolyl tRNA synthetase; Glutamyl tRNA synthetase; Glutamyl prolyl tRNA synthetase; Glutamyl tRNA synthetase; Glutamyl-tRNA synthetase; PARS; PIG 32; PIG32; Proliferation inducing gene 32 protein; Proliferation inducing protein 32; Proline tRNA ligase; Proline-tRNA ligase; Prolyl tRNA synthetase; Prolyl-tRNA synthetase; QARS; QPRS; SYEP_HUMAN;
Accession No.	Swiss-Prot#:P07814NCBI Gene ID:2058
Uniprot	P07814
GeneID	2058;
Calculated MW	171
Concentration	1mg/ml
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20 °C.Stable for 12 months from date of receipt

Application Details

WB: 1:500-1:2000

Images



Western blot analysis of Phospho-EPRS (Tyr377) in lysates of Jurkat?, using Phospho-EPRS (Tyr377) Antibody(13065).

Note: This product is for in vitro research use only