

PBOV1 Antibody

Catalog No: #40366



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description


Product Name	PBOV1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	WB IHC
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total PBOV1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide corresponding to residues near the N terminal of human prostate and breast cancer overexpressed 1
Target Name	PBOV1
Other Names	UC28; UROC28; dJ171N11.2
Accession No.	Swiss-Prot:Q9GZY1Gene Accssion:NP_067648
Uniprot	Q9GZY1
GeneID	59351;
SDS-PAGE MW	16KD
Concentration	1.4mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN3, 40% Glycerol.
Storage	Store at -20°C

Application Details


Western blotting: 1:500-1:2000

Immunohistochemistry:1:10-1:50

Images



Gel: 12%SDS-PAGE
Lysate: 40ug Mouse spleen tissuePrimary antibody: 1/200
dilution
Secondary antibody dilution: 1/8000
Exposure time: 40 seconds



Immunohistochemical analysis of paraffin-embedded Human breast cancer tissue using #40366 at dilution 1/20.

Background

This intronless gene encodes a protein of unknown function. Its expression is up-regulated in some types of cancer, including prostate, breast, and bladder cancer.

Note: This product is for in vitro research use only