

VEGFR1 (Phospho-Ser1048) Antibody

Catalog No: #13103



Package Size: #13103-1 50ul #13103-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	VEGFR1 (Phospho-Ser1048) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody is from purified rabbit serum by affinity purification via sequential chromatography on phospho- and non-phospho-peptide affinity columns.
Applications	WB
Species Reactivity	Hu, Ms, Rt
Specificity	Phospho-VEGFR1 (Phospho-Ser1048) Antibody detects endogenous levels of VEGFR1 only when phosphorylated at Ser1048
Immunogen Description	A synthesized peptide derived from human VEGFR1 around the phosphorylation site of Ser1048
Target Name	VEGFR1
Other Names	EC 2.7.10.1; FLT 1; FLT; Flt-1; FLT1; Fms like tyrosine kinase 1; Fms related tyrosine kinase 1; Fms related tyrosine kinase 1 (vascular endothelial growth factor/vascular permeability factor receptor); Fms related tyrosine kinase 1 vascular endothelial growth factor/vascular permeability factor receptor; Fms-like tyrosine kinase 1; FRT; Soluble VEGF receptor 1 14; Soluble VEGFR1 variant 2; Soluble VEGFR1 variant 21; Tyrosine protein kinase FRT; Tyrosine protein kinase receptor FLT; Tyrosine-protein kinase FRT; Tyrosine-protein kinase receptor FLT; Vascular endothelial growth factor receptor 1; Vascular endothelial growth factor vascular permeability factor receptor; Vascular permeability factor receptor 1; Vascular permeability factor receptor; VEGFR 1; VEGFR-1; VEGFR1; VGFR1_HUMAN;
Accession No.	Swiss-Prot#:P17948NCBI Gene ID:2321
Uniprot	P17948
GeneID	2321;
Calculated MW	151kda
Concentration	1mg/ml
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20 °C.Stable for 12 months from date of receipt

Application Details

WB: 1:500-1:2000

Images



Western blot analysis of Phospho-VEGFR1 (Ser1048) in lysates of K562 etoposide 25uM 24h, using Phospho-VEGFR1 (Ser1048) Antibody(13103).

Note: This product is for in vitro research use only