

BTK (Phospho-Tyr223/225) Antibody

Catalog No: #13202

Package Size: #13202-1 50ul #13202-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	BTK (Phospho-Tyr223/225) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody is from purified rabbit serum by affinity purification via sequential chromatography on phospho- and non-phospho-peptide affinity columns.
Applications	WB
Species Reactivity	Hu, Ms, Rt
Specificity	Phospho-BTK (Phospho-Tyr223/225) Antibody detects endogenous levels of BTK only when phosphorylated at Tyr223/225
Immunogen Description	A synthesized peptide derived from human BTK around the phosphorylation site of Tyr223/225
Target Name	BTK
Other Names	Agammaglobulinaemia tyrosine kinase; AGMX 1; AGMX1; AT; ATK; B cell progenitor kinase; B-cell progenitor kinase; BPK; Bruton agammaglobulinemia tyrosine kinase; Bruton tyrosine kinase; Brutons Tyrosine Kinase; Btk; BTK_HUMAN; dominant-negative kinase-deficient Brutons tyrosine kinase; IMD 1; IMD1; MGC126261; MGC126262; OTTHUMP00000063593; PSCTK 1; PSCTK1; truncated Bruton agammaglobulinemia tyrosine kinase; Tyrosine protein kinase BTK; Tyrosine-protein kinase BTK; tyrosine-protein kinase BTK isoform (lacking exon 14; XLA;
Accession No.	Swiss-Prot#:Q06187NCBI Gene ID:695
Uniprot	Q06187
GeneID	695;
Calculated MW	77
Concentration	1mg/ml
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20 °C.Stable for 12 months from date of receipt

Application Details

WB: 1:500-1:2000

Images



Western blot analysis of Phospho-BTK (Tyr223/225) in lysates of HeLa , using Phospho-BTK (Tyr223/225) Antibody(13202).

Note: This product is for in vitro research use only