

Ezh2 (Phospho-Ser363/366) Antibody

Catalog No: #13246



Package Size: #13246-1 50ul #13246-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Ezh2 (Phospho-Ser363/366) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody is from purified rabbit serum by affinity purification via sequential chromatography on phospho- and non-phospho-peptide affinity columns.
Applications	WB
Species Reactivity	Hu, Ms, Rt
Specificity	Phospho-Ezh2 (Phospho-Ser363/366) Antibody detects endogenous levels of Ezh2 only when phosphorylated at Ser363/366
Immunogen Type	peptide
Immunogen Description	A synthesized peptide derived from human Ezh2 around the phosphorylation site of Ser363/366
Target Name	Ezh2
Other Names	Enhancer of zeste 2; enhancer of zeste 2 polycomb repressive complex 2 subunit; Enhancer of zeste homolog 2 (Drosophila); Enhancer of zeste homolog 2; Enhancer of zeste, Drosophila, homolog 2; ENX 1; Enx 1h; ENX-1; ENX1; Enx1h; EZH 2; EZH1; EZH2; EZH2_HUMAN; EZH2b; Histone-lysine N-methyltransferase EZH2; KMT 6; KMT6; KMT6A; Lysine N-methyltransferase 6; MGC9169; WVS; WVS2;
Accession No.	Swiss-Prot#:Q15910NCBI Gene ID:2146
Uniprot	Q15910
GeneID	2146;
Calculated MW	85
Concentration	1mg/ml
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20 °C.Stable for 12 months from date of receipt

Application Details

WB: 1:500-1:2000

Images



Western blot analysis of Phospho-Ezh2 (Ser363/366) in lysates of Jurkat?, using Phospho-Ezh2 (Ser363/366) Antibody(13246).

Note: This product is for in vitro research use only