

Recombinant Human LR3 Insulin-Like Growth Factor I,LR3-IGF-1 (MG)

Catalog No: #AP74785

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74785-1 10ug #AP74785-2 50ug #AP74785-3 500ug #AP74785-4 1mg

Description

| | |
|-----------------|--|
| Product Name | Recombinant Human LR3 Insulin-Like Growth Factor I,LR3-IGF-1 (MG) |
| Host Species | E.coli |
| Purification | Greater than 95% as determined by reducing SDS-PAGE. |
| Specificity | Human |
| Accession No. | P05019 |
| Uniprot | P05019 |
| GeneID | 3479; |
| Calculated MW | 9.1kDa |
| Target Sequence | MFPAMPLSSLFVNGPRTLGAELVDALQFVCGDRGFYFNKPTGYGSSRRAPQTGIVDECCFRSCDLRRLEM YCAPLKPAKSA |
| Formulation | Lyophilized from a 0.2 µm filtered solution of 300mM HAc-NaAc, pH 6.5. |
| Storage | Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months. |

Note: This product is for in vitro research use only