

Recombinant Human Apoptosis-Inducing Factor 1 Mitochondrial,AIFM1 (N-6His)



Catalog No: #AP75004

Package Size: #AP75004-1 10ug #AP75004-2 50ug #AP75004-3 500ug #AP75004-4 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Apoptosis-Inducing Factor 1 Mitochondrial,AIFM1 (N-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O95831
Uniprot	O95831
GeneID	9131;
Calculated MW	56.2kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMEEVPQDKAPSHVPFLLIGGGTAFAAAARSIRARDPGARVLIVSEDPPELYM RPPLSKELWFSDDPNVTKTLRFKQWNGKERSIYFQPPSFYVSAQDLPHIENGGAVALTGKKVQLDVRDNMV KLNDGSQITYEKCLIATGGTPRSLAIDRAGAIEVKSRTTLFRKIGDFRSLKISREVKSITIIGGGFLGSELACALG RKARALGTEVIQLFPEKGNMGKILPEYLSNWTMEKVRREGVKVMPNAIVQSVGVSSGKLLIKLDGRKVVTDHI VAVGLEPNVELAKTGGLIDSDFGGFRVNAELQARSNIWVAGDAACFYDIKLGRRRVEHHDHAVVSGRLAG ENMTGAAKPYWHQSMFWSDLGPDVGYEAIGLVDSSLPTVGVFATAQDNPKSATEQSGTGIRSESETESE ASEITIPPSTPAVPQAPVQGEDYGGKGVIFYLRDKVVVGIVLWNIFNRMPIARKIIKDGEQHEDLNEVAKLFNIHED
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only