

Recombinant Human Multiple Inositol Polyphosphate Phosphatase 1,MINPP1 (C-6His)



Catalog No: #AP75223

Package Size: #AP75223-1 10ug #AP75223-2 50ug #AP75223-3 500ug #AP75223-4 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Multiple Inositol Polyphosphate Phosphatase 1,MINPP1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9UNW1
Uniprot	Q9UNW1
GeneID	9562;
Calculated MW	53.14kDa
Target Sequence	SLLEPRDPVASSLSPYFGTKTRYEDVNPVLLSGPEAPWRDPELLEGTCTPVQLVALIRHGTRYPTVKQIRKLRLQ LHGLLQARGSRDGGASSTGSRDLGAALADWPLWYADWMDGQLVEKGRQDMRQLALRLASLFPALFSRENY GRLRLITSSKHRCMDSSAAFLQGLWQHYHPGLPPPDVADMEFGPPTVNDKLMRFFDHCEKFLTEVEKNATAL YHVEAFKTGPEMQNILKKVAATLQVPVNDLNADLIQVAFFTCSFDLAIKGVKSPWCDVFDIDDAKVLEYLNDLK QYWKRGYGYTINSRSSCTLFQDIFQHLDKAVEQKQRSQPISSPVILQFGHAETLLPLLSLMGYFKDKEPLTAYN YKKQMHRKFRSGLIVPYASNLIFVLYHCENAKTPKEQFRVQMLLNEKVLPLAYSQETVSFYEDLKNHYKDILQS CQTSEECELARANSTSELDVDDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 10% Glycerol, pH 7.5.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only