



Recombinant Human High Mobility Group Protein B3, HMGB3 (C-6His)

Catalog No: #AP75401

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75401-1 10ug #AP75401-2 50ug #AP75401-3 500ug #AP75401-4 1mg

Description

Product Name	Recombinant Human High Mobility Group Protein B3, HMGB3 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O15347
Uniprot	O15347
GeneID	3149;
Calculated MW	24kDa
Target Sequence	MAKGDPKKPKGKMSAYAFFVQTCREEHKKKNPEVPVNFAEFSKKCSERWKTMSGKEKSKFDEMAKADKVR YDREMKDYGPAKGGKKKKDPNAPKRPPSGFFLCSEFRPKIKSTNPGISIGDVAKKLGEMWNNLNDSEKQPYI TKAAKLKEYEKDVADYKSKGKFDGAKGPAKVARKKVEEEDEEEEEEEEEEEEEDEVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only