

AGXT2L1 antibody

Catalog No: #22545

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	AGXT2L1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu Ms
Immunogen Type	Peptide-KLH
Immunogen Description	A synthesized peptide derived from human AGXT2L1, corresponding to a region within the internal amino acids.
Target Name	AGXT2L1
Accession No.	Swiss-Prot#: Q8TBG4NCBI Gene ID: 64850NCBI mRNA#: NM_031279NCBI Protein#: NP_112569
Uniprot	Q8TBG4
GeneID	64850;
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 56kd

Western blotting: 1:500-1:3000

Images



Western blot analysis of AGXT2L1 using HeLa whole cell lysates

Note: This product is for in vitro research use only