

Recombinant Human R-Spondin-3,RSPO3
(C-Fc-6His)

Catalog No: #AP75526



Package Size: #AP75526-1 10ug #AP75526-2 50ug #AP75526-3 500ug #AP75526-4 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human R-Spondin-3,RSPO3 (C-Fc-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9BXY4
Uniprot	Q9BXY4
GeneID	84870;
Calculated MW	56.3kDa
Target Sequence	QNASRGRQRMMHPNVSQGCQGGCATCSDYNGCLSCKPRLFFALERIGMKQIGVCLSSCPSGYYGTRYPDI NKCTKCKADCDCFCFNKFNCTKCKSGFYHLGLKCLDNCPEGLEANNHTMECVSIVHCEVSEWNPWSPCTKKG KTCGFKRGTTETRVREIIQHPSAKGNLCPPTNETRKCTVQRKKCQKGERGKKGRERKRKKPNKGESKEAIPDS KSLESSKEIPEQRENKQQQKKRKVQDKQKSVSVSTVHVDDIEGRMDEPKSCDKTHTCPPCPAPELLGGPSVF LFPPKPKDTLMISRTPEVTCVVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDW LNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQ PENNYKTTTPVLDSGDSFFLYSKLTVDKSRWQQGNVFCSSVMHEALHNHYTQKSLSLSPGKHHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only