

Recombinant Mouse Interleukin-6,IL-6

Catalog No: #AP75858



Package Size: #AP75858-1 10ug #AP75858-2 50ug #AP75858-3 500ug #AP75858-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse Interleukin-6,IL-6
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P08505
Uniprot	P08505
GeneID	16193;
Calculated MW	21.8kDa
Target Sequence	MFPTSQVRRGDFTEDTTPNRPVYTTSQVGGLITHVLWEIVEMRKELCNGNSDCMNNDDALAENNLKLPEIQR NDGCYQTGYNQEICLLKISSGLLEYHSYLEYMKNNLKDKNKKDKARVLQRDTETLIHIFNQEVKDLHKIVLPTPIS NALLTDKLESQKEWLRTKTIQFILKSLEEFKVTLRSTRQT
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only