

Recombinant Human IA2,PTPRN (Met687-Gln979, N-6His)

Catalog No: #AP75969

Package Size: #AP75969-1 10ug #AP75969-2 50ug #AP75969-3 500ug #AP75969-4 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human IA2,PTPRN (Met687-Gln979, N-6His)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q16849
Uniprot	Q16849
GeneID	5798;
Calculated MW	35.6kDa
Target Sequence	MGSSHHHHHHSSGLVPRGSHMDISTGHMILAYMEDHLRNRDRLAKEWQALCAYQAEPNTCATAQGEGNIKK NRHPDFLPYDHARIKLKVESSPSRSDYINASPIIEHDPRMPAYIATQGPLSHTIADFQMWVWESGCTVIVMLTPL VEDGVKQCDDRYWPDEGASLYHVYEVNLVSEHIWCEDFLVRSFYLNKVTQETRTLTFHFLSWPAEGTPAST RPLLDFFRRKVNKCYRGRSCPIIVHCSDGAGRTGTYLIDMVLNRMAKGVKEIDIAATLEHVRDQRPGLVRSKDQ FEFALTAVAAEEVNAILKALPQ
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris,150mM NaCl,pH8.0.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only