

Recombinant Human Thrombopoietin,TPO (N, C-6His)

Catalog No: #AP76129

Package Size: #AP76129-1 10ug #AP76129-2 50ug #AP76129-3 500ug #AP76129-4 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Thrombopoietin,TPO (N, C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P40225
Uniprot	P40225
GeneID	7066;
Calculated MW	37.3kDa
Target Sequence	HHHHHHSPAPPACDLRVLSKLLRDSHVLHSRLSQCEVHPLPTPVLLPAVDFSLGEWKTQMEETKAQDILGAV TLLLEGVMAARGQLGPTCLSSLLGQLSGQVRLLLGALQSLLGTQLPPQGRTTAHKDPNAIFLSFQHLLRGKVR FLMLVGGSTLCVRRAPPTTAVPSRTSLVLTNLNPNRTSGLLETNFTASARTTGSGLLKWQQGFRAKIPGLLN QTSRSLDQIPGYLNRIHELLNGTRGLFPGPSRRTLGAPODISSGTSDTGSLPPNLQPGYSPSPTHPPTGQYTLFP LPPTLPTPVVQLHPLLDPDSAPTPTPTSPLLNTSYTHSQNLSQEGVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM Tris,150mM NaCl,pH 8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only