

Recombinant Human Brain Natriuretic Peptide, BNP (N-6His-Flag)

Catalog No: #AP76236

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP76236-1 10ug #AP76236-2 50ug #AP76236-3 500ug #AP76236-4 1mg

Description

Product Name	Recombinant Human Brain Natriuretic Peptide, BNP (N-6His-Flag)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P16860
Uniprot	P16860
GeneID	4879;
Calculated MW	11kDa
Target Sequence	MNHNKVVHHHHHMDYKDDDDKHPLGSPGSASDLETSGLQEQRNHLQGKLSLELQVEQTSLEPLQESPRPTGV WKSREVATEGIRGHRKMLVLYTLRAPR
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only