

Recombinant Human Tryptase alpha,beta-1,TPSAB1
(C-6His)

Catalog No: #AP76354

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Package Size: #AP76354-1 10ug #AP76354-2 50ug #AP76354-3 500ug #AP76354-4 1mg

Description

Product Name	Recombinant Human Tryptase alpha,beta-1,TPSAB1 (C-6His)
Host Species	Human Cells
Purification	Greater than 90% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q15661
Uniprot	Q15661
GeneID	7177;
Calculated MW	28.8kDa
Target Sequence	IVGGQEAPRSKWPWQVSLRVRDRYWMHFCGGSLIHPQWVLTAAHCLGPDVKDLATLRVQLREQHLYYQDQ LLPVSRIIVHPQFYIIQTGADIALLEEPVNISSRVHTVMLPPASETFPPGMPWCWVTGWGDVDNDEPLPPFPPL KQVKVPIMENHICDAKYHLGAYTGDDVRIIRDDMLCAGNTRRDSCQGDSGGPLVCKVNGTWLQAGVVSWE GCAQPNRPGIYTRVTYYLDWIHHYVPKKPVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only