

Recombinant Human GTP-binding protein Di-Ras3(DIRAS3)

Catalog No: #AP76790

Package Size: #AP76790-1 20ug #AP76790-2 100ug #AP76790-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human GTP-binding protein Di-Ras3(DIRAS3)
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:1-229aaSequence Info:Full Length
Other Names	Distinct subgroup of the Ras family member 3 Rho-related GTP-binding protein Rho1
Accession No.	O95661
Uniprot	O95661
GeneID	9077;
Calculated MW	52.5 kDa
Tag Info	N-terminal GST-tagged
Target Sequence	MGNASFGSKEQKLLKRLRLLPALLILRAFKPHRKIRDYRVVVVGTAGVGKSTLLHKWASGNFRHEYLP TIENTY CQLLGCSHGVL SLHITDSKSGDGNRALQRHVIARGHAFVLVYSVTKKETLEELKAFYELICKIKGNNLHKFPVL VGNKSDDTHREVALNDGATCAMEWNCAFMESISAKTDVNVQELFHMLLNYKKKPTTGLQEPEKKSQMPNTTE KLLDKC
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

References

"NOEY2 (ARHI), an imprinted putative tumor suppressor gene in ovarian and breast carcinomas."

Yu Y.H., Xu F.J., Peng H., Fang X., Zhao S., Li Y., Cuevas B., Kuo W.-L., Gray J.W., Siciliano M., Mills G.B., Bast R.C. Jr.

Proc. Natl. Acad. Sci. U.S.A. 96:214-219(1999) Research Topic:Others

Note: This product is for in vitro research use only