

# AmCyan Monoclonal Antibody

Catalog No: #T614



Package Size: #T614-1 50ul #T614-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	AmCyan Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	WB
Specificity	The antibody detects AmCyan & AmCyan Tag protein.
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western blotting: 1:5000

## Images



Western blot analysis of recombinant AmCyan protein, diluted at 1:5000.

Note: This product is for in vitro research use only