

EBFP Monoclonal Antibody

Catalog No: #T615



Package Size: #T615-1 50ul #T615-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

| | |
|---------------|--|
| Product Name | EBFP Monoclonal Antibody |
| Host Species | Mouse |
| Clonality | Monoclonal |
| Purification | The antibody was affinity-purified from antiserum by affinity-chromatography |
| Applications | WB |
| Specificity | The antibody detects EBFP & EBFP Tag protein. |
| Concentration | 1mg/ml |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Store at -20°C/1 year |

Application Details

Western blotting: 1:5000

Images



Western blot analysis of recombinant EBFP protein, diluted at 1:5000.

Note: This product is for in vitro research use only