

FH Monoclonal Antibody

Catalog No: #40442

Package Size: #40442-1 50ul #40442-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	FH Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	WB IF
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous FH proteins.
Target Name	FH
Accession No.	Swiss-Prot: P07954NCBI Gene ID: 2271
Uniprot	P07954
GeneID	2271;
SDS-PAGE MW	50kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:3000

Images



Western blot analysis of 1) 293T, 2) HepG2, 3) HeLa, diluted at 1:3000.

Note: This product is for in vitro research use only