

# Rabbit IgG (Fc specific, Heavy Chain) Monoclonal Antibody

Catalog No: #40497

Package Size: #40497-1 50ul #40497-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rabbit IgG (Fc specific, Heavy Chain) Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	WB
Species Reactivity	Rabbit
Specificity	The antibody reacts with the Fc portion of rabbit IgG heavy chain but not with the Fab portion of rabbit immunoglobulins. The antibody may cross-react with immunoglobulins from other species.
SDS-PAGE MW	55kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western blotting: 1:2000-5000

## Images



Western blot analysis of Reduced Rabbit IgG, H+L, diluted at 1:10000.

Note: This product is for in vitro research use only