

Cleaved-MMP-10 (F99) Polyclonal Antibody

Catalog No: #40514

Package Size: #40514-1 50ul #40514-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cleaved-MMP-10 (F99) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	Cleaved-MMP-10 (F99) Polyclonal Antibody detects endogenous levels of fragment of activated MMP-10 protein resulting from cleavage adjacent to F99.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human MMP-10.
Target Name	Cleaved-MMP-10
Other Names	MMP10; STMY2; Stromelysin-2; SL-2; Matrix metalloproteinase-10; MMP-10; Transin-2
Accession No.	Swiss-Prot: P09238NCBI Gene ID: 4319
Uniprot	P09238
GeneID	4319;
SDS-PAGE MW	54kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.


ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of HepG2 cells using Cleaved-MMP-10 (F99) Polyclonal Antibody



Western Blot analysis of NIH-3T3 cells using
Cleaved-MMP-10 (F99) Polyclonal Antibody

Note: This product is for in vitro research use only