

ABCC13 Polyclonal Antibody

Catalog No: #40534

Package Size: #40534-1 50ul #40534-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ABCC13 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu
Specificity	ABCC13 Polyclonal Antibody detects endogenous levels of ABCC13 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human ABCC13.
Target Name	ABCC13
Other Names	ABCC13; C21orf73; PRED6; Putative ATP-binding cassette sub-family C member 13
Accession No.	Swiss-Prot: Q9NSE7NCBI Gene ID: 150000
Uniprot	Q9NSE7
GeneID	150000;
SDS-PAGE MW	31kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of HT29 cells using ABCC13 Polyclonal Antibody



Western Blot analysis of HELA cells using ABCC13 Polyclonal Antibody

Note: This product is for in vitro research use only