

ACE1 Polyclonal Antibody

Catalog No: #40540

Package Size: #40540-1 50ul #40540-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ACE1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	ACE1 Polyclonal Antibody detects endogenous levels of ACE1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human ACE1.
Target Name	ACE1
Other Names	ACE; DCP; DCP1; Angiotensin-converting enzyme; ACE; Dipeptidyl carboxypeptidase I; Kininase II; CD antigen CD143
Accession No.	Swiss-Prot: P12821NCBI Gene ID: 1636
Uniprot	P12821
GeneID	1636;
SDS-PAGE MW	165kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of HuvEc cells using ACE1 Polyclonal Antibody

Note: This product is for in vitro research use only