

AChR α 10 Polyclonal Antibody

Catalog No: #40542

Package Size: #40542-1 50ul #40542-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description


Product Name	AChR α 10 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	AChR α 10 Polyclonal Antibody detects endogenous levels of AChR α 10 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human AChR α 10.
Target Name	AChR α 10
Other Names	CHRNA10; NACHRA10; Neuronal acetylcholine receptor subunit alpha-10; Nicotinic acetylcholine receptor subunit alpha-10; NACHR alpha-10
Accession No.	Swiss-Prot: Q9GZZ6NCBI Gene ID: 57053
Uniprot	Q9GZZ6
GeneID	57053;
SDS-PAGE MW	50kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Images

Western Blot analysis of A549 cells using AChR α 10 Polyclonal Antibody



Western Blot analysis of HuvEc MCF7 cells using AChRα10 Polyclonal Antibody

Note: This product is for in vitro research use only