

Aldolase A Polyclonal Antibody

Catalog No: #40581

Package Size: #40581-1 50ul #40581-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Aldolase A Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	Aldolase A Polyclonal Antibody detects endogenous levels of Aldolase A protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human Aldolase A.
Target Name	Aldolase A
Other Names	ALDOA; ALDA; Fructose-bisphosphate aldolase A; Lung cancer antigen NY-LU-1; Muscle-type aldolase
Accession No.	Swiss-Prot: P04075NCBI Gene ID: 226
Uniprot	P04075
GeneID	226;
SDS-PAGE MW	39kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunofluorescence: 1/200 - 1/1000.


ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of HT29 cells using Aldolase A Polyclonal Antibody



Western Blot analysis of NIH-3T3 K562 cells using Aldolase A Polyclonal Antibody

Note: This product is for in vitro research use only