

Annexin III Polyclonal Antibody

Catalog No: #40593



Package Size: #40593-1 50ul #40593-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Annexin III Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms
Specificity	Annexin III Polyclonal Antibody detects endogenous levels of Annexin III protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Annexin III.
Target Name	Annexin III
Other Names	ANXA3; ANX3; Annexin A3; 35-alpha calcimedin; Annexin III; Annexin-3; Inositol 1; 2-cyclic phosphate 2-phosphohydrolase; Lipocortin III; Placental anticoagulant protein III; PAP-III
Accession No.	Swiss-Prot: P12429NCBI Gene ID: 306
Uniprot	P12429
GeneID	306;
SDS-PAGE MW	36kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


ELISA: 1/5000.

Not yet tested in other applications.

Images



Western Blot analysis of HepG2 cells using Annexin III Polyclonal Antibody



Western Blot analysis of K562 cells using Annexin III Polyclonal Antibody

Note: This product is for in vitro research use only