

Annexin VI Polyclonal Antibody

Catalog No: #40594

Package Size: #40594-1 50ul #40594-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Annexin VI Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu
Specificity	Annexin VI Polyclonal Antibody detects endogenous levels of Annexin VI protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human Annexin VI.
Target Name	Annexin VI
Other Names	ANXA6; ANX6; Annexin A6; 67 kDa calelectrin; Annexin VI; Annexin-6; Calphobindin-II; CPB-II; Chromobindin-20; Lipocortin VI; Protein III; p68; p70
Accession No.	Swiss-Prot: P08133NCBI Gene ID: 309
Uniprot	P08133
GeneID	309;
SDS-PAGE MW	75kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.


ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of 293 cells using Annexin VI Polyclonal Antibody



Western Blot analysis of 293 cells using Annexin VI Polyclonal Antibody

Note: This product is for in vitro research use only