

Apaf-1-ALT Polyclonal Antibody

Catalog No: #40600

Package Size: #40600-1 50ul #40600-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Apaf-1-ALT Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu
Specificity	Apaf-1-ALT Polyclonal Antibody detects endogenous levels of Apaf-1-ALT protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Apaf-1-ALT.
Target Name	Apaf-1-ALT
Other Names	APAF1-interacting protein; APIP; APIP2; CGI-29 protein; MMRP19; dJ179L10.2; likely ortholog of mouse monocyte macrophage 19
Accession No.	Swiss-Prot: O14727-6NCBI Gene ID: 317
Uniprot	O14727
GeneID	317;
SDS-PAGE MW	27kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.


ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of COLO205 cells using Apaf-1-ALT Polyclonal Antibody



Western Blot analysis of Jurkat cells using Apaf-1-ALT Polyclonal Antibody

Note: This product is for in vitro research use only