

APC Polyclonal Antibody

Catalog No: #40601

Package Size: #40601-1 50ul #40601-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	APC Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	APC Polyclonal Antibody detects endogenous levels of APC protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human APC.
Target Name	APC
Other Names	APC; DP2.5; Adenomatous polyposis coli protein; Protein APC; Deleted in polyposis 2.5
Accession No.	Swiss-Prot: P25054NCBI Gene ID: 324
Uniprot	P25054
GeneID	324;
SDS-PAGE MW	310kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.


ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of COLO205 cells using APC Polyclonal Antibody



Western Blot analysis of RAT-MUSCLE cells using APC Polyclonal Antibody

Note: This product is for in vitro research use only