

AT1 Polyclonal Antibody

Catalog No: #40620



Package Size: #40620-1 50ul #40620-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	AT1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	AT1 Polyclonal Antibody detects endogenous levels of AT1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human AT1.
Target Name	AT1
Other Names	AGTR1; AGTR1A; AGTR1B; AT2R1; AT2R1B; Type-1 angiotensin II receptor; AT1AR; AT1BR; Angiotensin II type-1 receptor; AT1
Accession No.	Swiss-Prot: P30556NCBI Gene ID: 185
Uniprot	P30556
GeneID	185;
SDS-PAGE MW	41kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of HEPG2 cells using AT1 Polyclonal Antibody

Note: This product is for in vitro research use only