

BAF250b Polyclonal Antibody

Catalog No: #40632



Package Size: #40632-1 50ul #40632-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	BAF250b Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Rt
Specificity	BAF250b Polyclonal Antibody detects endogenous levels of BAF250b protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human BAF250b.
Target Name	BAF250b
Other Names	ARID1B; BAF250B; DAN15; KIAA1235; OSA2; AT-rich interactive domain-containing protein 1B; ARID domain-containing protein 1B; BRG1-associated factor 250b; BAF250B; BRG1-binding protein hELD/OSA1; Osa homolog 2; hOsa2; p250R
Accession No.	Swiss-Prot: Q8NFD5NCBI Gene ID: 57492
Uniprot	Q8NFD5
GeneID	57492;
SDS-PAGE MW	170kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of 3T3 cells using BAF250b Polyclonal Antibody

Note: This product is for in vitro research use only