

Bcl-6 Polyclonal Antibody

Catalog No: #40640



Package Size: #40640-1 50ul #40640-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bcl-6 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	Bcl-6 Polyclonal Antibody detects endogenous levels of Bcl-6 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Bcl-6 around the non-phosphorylation site of S333.
Target Name	Bcl-6
Other Names	BCL6; BCL5; LAZ3; ZBTB27; ZNF51; B-cell lymphoma 6 protein; BCL-6; B-cell lymphoma 5 protein; BCL-5; Protein LAZ-3; Zinc finger and BTB domain-containing protein 27; Zinc finger protein 51
Accession No.	Swiss-Prot: P41182NCBI Gene ID: 604
Uniprot	P41182
GeneID	604;
SDS-PAGE MW	90kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of 293 cells using Bcl-6 Polyclonal Antibody



Western Blot analysis of mouse-heart NIH-3T3 cells using Bcl-6 Polyclonal Antibody

Note: This product is for in vitro research use only