

Caspase-8 Polyclonal Antibody

Catalog No: #40674



Package Size: #40674-1 50ul #40674-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Caspase-8 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	Caspase-8 Polyclonal Antibody detects endogenous levels of Caspase-8 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Caspase-8 around the non-phosphorylation site of S347.
Target Name	Caspase-8
Other Names	CASP8; MCH5; Caspase-8; CASP-8; Apoptotic cysteine protease; Apoptotic protease Mch-5; CAP4; FADD-homologous ICE/ced-3-like protease; FADD-like ICE; FLICE; ICE-like apoptotic protease 5; MORT1-associated ced-3 homolog; MACH
Accession No.	Swiss-Prot: Q14790NCBI Gene ID: 841
Uniprot	Q14790
GeneID	841;
SDS-PAGE MW	55kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


ELISA: 1/5000.

Not yet tested in other applications.

Images



Western Blot analysis of 3T3 cells using Caspase-8 Polyclonal Antibody



Western Blot analysis of LOVO cells using Caspase-8 Polyclonal Antibody

Note: This product is for in vitro research use only