

## CBP Polyclonal Antibody

Catalog No: #40686



Package Size: #40686-1 50ul #40686-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                                                                                       |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name          | CBP Polyclonal Antibody                                                                                               |
| Host Species          | Rabbit                                                                                                                |
| Clonality             | Polyclonal                                                                                                            |
| Purification          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Applications          | WB IHC IF ELISA                                                                                                       |
| Species Reactivity    | Hu Ms Rt                                                                                                              |
| Specificity           | CBP Polyclonal Antibody detects endogenous levels of CBP protein.                                                     |
| Immunogen Type        | peptide                                                                                                               |
| Immunogen Description | Synthesized peptide derived from human CBP around the non-acetylation site of K1535.                                  |
| Target Name           | CBP                                                                                                                   |
| Other Names           | CREBBP; CBP; CREB-binding protein                                                                                     |
| Accession No.         | Swiss-Prot: Q92793NCBI Gene ID: 1387                                                                                  |
| Uniprot               | Q92793                                                                                                                |
| GeneID                | 1387;                                                                                                                 |
| SDS-PAGE MW           | 265kd                                                                                                                 |
| Concentration         | 1mg/ml                                                                                                                |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                               |
| Storage               | Store at -20°C/1 year                                                                                                 |

## Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.


ELISA: 1/20000.

Not yet tested in other applications.

## Images



Western Blot analysis of HT29 cells using CBP Polyclonal Antibody



Western Blot analysis of AD293 HELA cells using CBP Polyclonal Antibody

Note: This product is for in vitro research use only