

CD83 Polyclonal Antibody

Catalog No: #40707

Package Size: #40707-1 50ul #40707-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD83 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Species Reactivity	Hu Ms Rt
Specificity	CD83 Polyclonal Antibody detects endogenous levels of CD83 protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human CD83.
Target Name	CD83
Other Names	CD83; CD83 antigen; hCD83; B-cell activation protein; Cell surface protein HB15; CD antigen CD83
Accession No.	Swiss-Prot: Q01151NCBI Gene ID: 9308
Uniprot	Q01151
GeneID	9308;
SDS-PAGE MW	23kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of Jurkat cells using CD83 Polyclonal Antibody

Note: This product is for in vitro research use only