

CD9 Polyclonal Antibody

Catalog No: #40708

Package Size: #40708-1 50ul #40708-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD9 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	CD9 Polyclonal Antibody detects endogenous levels of CD9 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human CD9.
Target Name	CD9
Other Names	CD9; MIC3; TSPAN29; GIG2; CD9 antigen; 5H9 antigen; Cell growth-inhibiting gene 2 protein; Leukocyte antigen MIC3; Motility-related protein; MRP-1; Tetraspanin-29; Tspan-29; p24; CD antigen CD9
Accession No.	Swiss-Prot: P21926NCBI Gene ID: 928
Uniprot	P21926
GeneID	928;
SDS-PAGE MW	25kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.


ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of 3T3 cells using CD9 Polyclonal Antibody



Western Blot analysis of L929 cells using CD9 Polyclonal Antibody

Note: This product is for in vitro research use only