

CDKN3 Polyclonal Antibody

Catalog No: #40724

Package Size: #40724-1 50ul #40724-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	CDKN3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC
Species Reactivity	Hu Ms
Specificity	CDKN3 Polyclonal Antibody detects endogenous levels of CDKN3 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human CDKN3.
Target Name	CDKN3
Other Names	CDKN3; CDI1; CIP2; KAP; Cyclin-dependent kinase inhibitor 3; CDK2-associated dual-specificity phosphatase; Cyclin-dependent kinase interactor 1; Cyclin-dependent kinase-interacting protein 2; Kinase-associated phosphatase
Accession No.	Swiss-Prot: Q16667NCBI Gene ID: 1033
Uniprot	Q16667
GeneID	1033;
SDS-PAGE MW	23kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

Images



Western Blot analysis of SH-SY5Y cells using CDKN3 Polyclonal Antibody

Note: This product is for in vitro research use only